



# *In vitro* gene editing - RNA purification

## Purification of sgRNA with Quick-RNA™ Microprep Kit - 10 min

### Purification of sgRNA - Needed Materials:

- Quick-RNA™ Microprep Kit (<https://www.zymoresearch.com/products/quick-rna-microprep-kit>);
- Nuclease-free H<sub>2</sub>O (i.e., DNase- and RNase- free water);
- 100% ethanol;
- RNase-free tubes and filter pipette tips;
- Micropipettes.

### Purification of sgRNA - Protocol:

- a) Add an equal volume of **RNA Lysis Buffer** to the sample (50µL of buffer per 50µL of sample) and mix.
- b) Add two volumes of 100% ethanol (200µL of ethanol per 50µL of initial sample) and mix.



- c) Transfer the mixture to the **Spin Column** in a **Collection Tube** and centrifuge for 1 min at 13 000 x g. Discard flow-through;
- d) Add 400µL of **RNA Prep Buffer** to the column and centrifuge for 1 min at 13 000 x g. Discard flow-through.
- e) Add 700µL of **RNA Wash Buffer** to the column and centrifuge for 1 min at 13 000 x g. Discard flow-through.
- f) Add 400µL of **RNA Wash Buffer** to the column and centrifuge at 13 000 x g for 2 min to ensure complete removal of the wash buffer. Carefully, transfer the column into a nuclease-free tube.
- g) Add 20µL of Nuclease-free H<sub>2</sub>O directly to the column matrix and centrifuge for 1 min at 13 000 x g. Use immediately, or store RNA at -80°C.
- h) Analyze purified sgRNA by gel electrophoresis.

