

In vitro cleavage - 3h

In vitro cleavage of DNA using CRISPR/Cas9 - 3.00h

The objective of this step is to make the CRISPR/Cas9 reaction, run the *in vitro* digestion in a gel, interpret the results and discuss them.

SUMMARY OF SEQUENCES

- sgHL3 target site and sgRNA sequence
- sgHL4 target site and sgRNA sequence

In vitro cleavage of DNA using CRISPR/Cas9 - Needed Materials:

- Nuclease-free water;
- 10x Reaction Buffer;
- Synthetized sgRNA (sgHL3 or sgHL4);
- 500 ng/ μ L Cas9 protein;
- 100 ng/ μ L Target DNA (698bp fragment containing the target site);
- 20 mg/ml Proteinase K;
- RNase-free tubes and filter pipette tips;
- Micropipettes;
- 37°C / 65°C heat block/incubator.



***In vitro* cleavage of DNA using CRISPR/Cas9 - Protocol:**


- a)** Perform the DNA digestion reaction by assembling the following components:

| Reagent | Volume (μL) |
|--------------------------|--------------------|
| Nuclease-free water | 5 |
| 10x BSA (10% BSA) | 1.5 |
| 10x Reaction Buffer | 1.5 |
| sgRNA | 5 |
| Cas9 protein (500 ng/μL) | 1 |
| Target DNA (100 ng/μL) | 1 |
| Final Volume (μL) | 15 |

- b)** Mix gently and incubate at 37°C for 60-90 min.
- c)** Add 0.5μL of Proteinase K (20 mg/ml) and incubate at 65°C for 15 min to stop the reaction and release the DNA from the Cas9 endonuclease.
Proteinase K degrades proteins and will destroy the Cas9 protein, thereby releasing the bound DNA and RNA fragments.
- d)** Add an appropriate volume of Loading Dye and run the samples on a 1% agarose gel alongside 3μL of GeneRuler 1 kb DNA Ladder (refer to step [V. Agarose Gel Electrophoresis](#) for full protocol).

NOTE: The DNA fragments will not migrate in the gel unless released from the Cas9 protein during the previous step.



**Expected Results:**

- 100% cleavage efficiency: two bands from cleaved template
- <100% cleavage efficiency: three bands (two from cleaved template and one uncut)
- 0% cleavage efficiency: single band from uncut template.



SUMMARY OF SEQUENCES:

sgHL3 target site and sgRNA sequence

DESIGN OF sgRNA (sgHL3)

Target DNA sequence (698 bp):

5' -

CTTCTCCTCGACGTGCACGGACTGGACCTCGCCGAGCAGGATGATGTGGTCGCCCCGCCGGGTAGCGCTCG
TGGACCGTGCACCTCGACGACCGCGACCGCTCCGTCGAGCACGGTCGCTCCCCGCGCGGTACGGACGAAC
TCCCCGCCCCGGAACCTGTCCGCGGACTTGC GCGCGAAGCGCATGGCCAGGTCCGTTGTGGTCCCTCGCGCA
GCACGCTCACCGCGAACTCGCCGCAACTGTGCAACACCGGGAAGGAGTTGGCCGTACGAGCCAGGCAGA
CCAGTGCCAGTGGCGGCTCCATCGAGACGGACACGAACGAACTGGCGGTGAAACCGTGC GGGACTCCCC
CGCGGTATGGGCGGTGACGAGCGCCACCCCGGCCGACCCGGGCCATGGC**GTCCCGGAGCATTCCCT**
GGTCGGCTGCCATCTTGAACCTCCCTAGGCGAGGCAGGTGGGCAGCTGCCACCGTTGTAGGCGACCAGC
TCCGCGAAGATCGCCATGGCGTCACGCACGTGCCCGCCAGCAGCTCCTCGTGGGCCCGGTGCAGGTTCCG
CGACGAGGCGCTCCGCGTCCAGCAGCCCGCGAAGGCGCCGAGGTCCGCCCGGCGCGGCGCCTGCAAC
GGTGTAGCCGCGGTGACCGCATCGGCGGCGAGGCGCTGGACCCAGCGCAGATAGTCCCGGTTGGCG
TCGAT

- 3'

1. Select target site (use of a target DNA selection program is recommended)

sgRNA (sgHL3) target site: 5' - **GTCCCGGAGCATTCCCTGGT** - 3'

(At least one "**G**" is necessary for RNA transcription; if no "G" is present at the 5' end of the target site sequence, add one "G")

2. Append T7 RNA Polymerase promoter sequence to the 5' end:

5' - TTCTAATACGACTCACTATA**GTCCCGGAGCATTCCCTGGT** - 3'

3. Append Scaffold overlap sequence to the 3' end:

5' - TTCTAATACGACTCACTATA**GTCCCGGAGCATTCCCTGGT**GTTTTAGAGCTAGA - 3'

7 RNA polymerase

Target site

Scaffold overlap sequence



SYNTHESIS OF sgRNA (sgHL3)

sgHL3 oligo:

5' - TTCTAATACGACTCACTATAGTCCCGGAGCATTCCCTGGTGTTTTAGAGCTAGA - 3'

Scaffold oligo:

5' - AAAAGCACCGACTCGGTGCCACTTTTTCAAGTTGATAACGGACTAGCCTTATTTAA
CTTGCTATTTCTAGCTCTAAAAC - 3'*The 14 nt sequence at the 3' end (purple underlined) is complementary to the Scaffold overlap sequence at the 3' end of the sgRNA oligo (purple)***1. Generate sgHL3 dsDNA template (120 bp):***Amplicon sequence that results from PCR amplification using sgHL3 oligo + Scaffold oligo (no dsDNA template is needed):*5' - TTCTAATACGACTCACTATAGTCCCGGAGCATTCCCTGGTGTTTTAGAGCTAGAAA
TAGCAAGTTAAAATAAGGCTAGTCCGTTATCAACTTGAAAAAGTGGCACCGAGTCGGTGCTTTT - 3'**2. Transcribe RNA from dsDNA template using T7 RNA polymerase***The T7 RNA polymerase starts transcription at the underlined **G** in the promoter sequence (see above) and transcribes RNA using the opposite strand as a template from 5'→3', i.e. The first base in the transcript will be G.***sgHL3 (100 bp):**5' - GUCCCGGAGCAU**UCCUGGU**GUUUUAGAGCUAGAAAUAGCAAGUUAAAAUAAGG
CUAGUCCGUUAUCAACUUGAAAAAGUGGCACCGAGUCGGUGCUUUU - 3'*The resulting sgRNA contains the desired 20 nt sequence (orange) complementary to the target DNA sequence.*

IN VITRO CRISPR/CAS9 (sgHL3)**Target DNA sequence (698 bp):**

5' -

CTTCTCCTCGACGTGCACGGACTGGACCTCGCCGAGCAGGATGATGTGGTCGCCCCGCCGGGTAGCGCTCG
 TGGACCGTGCACACTCGACGACCGCGACCGCTCCGTTCGAGCACGGTCGCTCCCCGCGCGGTACGGACGAAC
 TCCCCGCCCGCGAACTTGTCCGCGGACTTGCGCGCGAAGCGCATGGCCAGGTCCGTGTGGTCCTCGCGCA
 GCACGCTCACCGCGAACTCGCCGCAACTGTGCAACACCGGGAAGGAGTTGGCCGTACGAGCCAGGCAGA
 CCAGTGCCAGTGGCGGCTCCATCGAGACGGACACGAACGAACTGGCGGTGAAACCGTGCGGGACTCCCC
 CGCGGTTCATGGGCGGTGACGAGCGCCACCCCGGCCGACCCGGGCCATGGC**GTCCCGGAGCATTCCCT**
GGTCGGCTGCCATCTTCGAACTCCCTAGGCGAGGCAGGTGGGCAGCTGCCACCGTTGTAGG
 CGACCAGCTCCGCGAAGATCGCCATGGCGTCACGCACGTGCCCGCCAGCAGCTCCTCGTGGGCCCGGT
 GCAGGTTTCGCGACGAGGCGCTCCGCGTCCAGCAGCCCGGCGAAGGCGCCGAGTCCGCCCGGCGCGCG
 GCCTGCAACGGTGTGACCCGGCGGTGACCCGCATCGGCGGCGAGGCGCTGGACCCAGCGCAGATAGTCC
 CGGTTGGCGTCGAT

- 3'

Cas9 protein complexed with sgRNA binds to complementary DNA sequence on the target strand after recognition of the PAM site (5'-NGG-3', in red). The Cas9 nuclease cuts 3 nt upstream of the PAM site (red line).

Cleaved sequence (416 + 282 bp):

5' -

CTTCTCCTCGACGTGCACGGACTGGACCTCGCCGAGCAGGATGATGTGGTCGCCCCGCCGGGTAGCGCTCG
 TGGACCGTGCACACTCGACGACCGCGACCGCTCCGTTCGAGCACGGTCGCTCCCCGCGCGGTACGGACGAAC
 TCCCCGCCCGCGAACTTGTCCGCGGACTTGCGCGCGAAGCGCATGGCCAGGTCCGTGTGGTCCTCGCGCA
 GCACGCTCACCGCGAACTCGCCGCAACTGTGCAACACCGGGAAGGAGTTGGCCGTACGAGCCAGGCAGA
 CCAGTGCCAGTGGCGGCTCCATCGAGACGGACACGAACGAACTGGCGGTGAAACCGTGCGGGACTCCCC
 CGCGGTTCATGGGCGGTGACGAGCGCCACCCCGGCCGACCCGGGCCATGGC**GTCCCGGAGCATTCCCT**


- 3'

+

5' -

GGTCGGCTGCCATCTTCGAACTCCCTAGGCGAGGCAGGTGGGCAGCTGCCACCGTTGTAGGCGACCAGC
 TCCGCGAAGATCGCCATGGCGTCACGCACGTGCCCGCCAGCAGCTCCTCGTGGGCCCGGTGCAGGTTCC
 CGACGAGGCGCTCCGCGTCCAGCAGCCCGGCGAAGGCGCCGAGTCCGCCCGGCGCGCGGCCTGCAAC
 GGTGTCAGCCGGCGGTGACCCGCATCGGCGGCGAGGCGCTGGACCCAGCGCAGATAGTCCCGGTTGGCG
 TCGAT - 3'





SUMMARY OF SEQUENCES:

sgHL4 target site and sgRNA sequence

DESIGN OF sgRNA (sgHL4)

Target DNA sequence (698 bp):

5' -

CTTCTCCTCGACGTGCACGGACTGGACCTCGCCGAGCAGGATGATGTGGTCGCCCGCCGGGTAGCGCTCG
TGGACCGTGCACCTCGACGACCGCGACCGCTCCGTCGAGCACGGTCGCTCCCCGCGCGGTACGGACGAAC
TCCCCGCCCGCGAACTTGTCCGCGGACTTGC GCGCGAAGCGCATGGCCAGGTCCGTCGTGGTCCCTCGCGCA
GCACGCTCACCGCGAACTCGCCGCAACTGTGCAACACCGGGAAGGAGTTGGCCGTACGAG**GCCAGGCAGA**
CCAGTGCCAGTGGCGGCTCCATCGAGACGGACACGAACGAACTGGCGGTGAAACCGTGC GGGACTCCCC
CGCGGTCATGGGCGGTGACGAGCGCCACCCCGGCCGACCCGGGCCATGGCGTCCCGGAGCATTCCCT
GGTCGGCTGCCATCTTGAACCTCCCTAGGCGAGGCAGGTGGCAGCTGCCACCGTTGTAGGCGACCAGC
TCCGCGAAGATCGCCATGGCGTCACGCACGTGCCCGCCAGCAGCTCCTCGTGGGCCCGGTGCAGGTTCC
CGACGAGGCGCTCCGCGTCCAGCAGCCCGCGAAGGCGCCGAGGTCCGCCCGGCGCGGCGCCTGCAAC
GGTGTGACCCGCGGTCGACCGCATCGGCGGCGAGGCGCTGGACCCAGCGCAGATAGTCCCGGTTGGCG
TCGAT

- 3'

1. Select target site (use of a target DNA selection program is recommended)

sgRNA (sgHL4) target site: 5' - **GCCAGGCAGACCAGTGCCAG** - 3'

(At least one "**G**" is necessary for RNA transcription; if no "G" is present at the 5' end of the target site sequence, add one "G")

2. Append T7 RNA Polymerase promoter sequence to the 5' end:

5' - TTCTAATACGACTCACTATA**GCCAGGCAGACCAGTGCCAG** - 3'

3. Append Scaffold overlap sequence to the 3' end:

5' - TTCTAATACGACTCACTATA**GCCAGGCAGACCAGTGCCAG**GTTTTAGAGCTAGA - 3'

7 RNA polymerase

Target site

Scaffold overlap sequence



SYNTHESIS OF sgRNA (sgHL4)

sgHL4 oligo:

5' - TTCTAATACGACTCACTATAGCCAGGCAGACCAGTGCCAGGTTTTAGAGCTAGA - 3'

Scaffold oligo:

5' - AAAAGCACCGACTCGGTGCCACTTTTTCAAGTTGATAACGGACTAGCCTTATTTTAA
CTTGCTATTTCTAGCTCTAAAAC - 3'*The 14 nt sequence at the 3' end (purple underlined) is complementary to the Scaffold overlap sequence at the 3' end of the sgRNA oligo (purple)***1. Generate sgHL4 dsDNA template (120 bp):***Amplicon sequence that results from PCR amplification using sgRNA oligo + Scaffold oligo (no dsDNA template is needed):*5' - TTCTAATACGACTCACTATAGCCAGGCAGACCAGTGCCAGGTTTTAGAGCTAGAAA
TAGCAAGTTAAAATAAGGCTAGTCCGTTATCAACTTGAAAAAGTGGCACCGAGTCGGTGCTTTT - 3'**2. Transcribe RNA from dsDNA template using T7 RNA polymerase***The T7 RNA polymerase starts transcription at the underlined **G** in the promoter sequence (see above) and transcribes RNA using the opposite strand as a template from 5'→3', i.e. The first base in the transcript will be G.***sgHL3 (100 bp):**5' - GCCAGGCAGACCAGUGCCAGGUUUUAGAGCUAGAAAUAGCAAGUUAAAAUAAGG
CUAGUCCGUUAUCAACUUGAAAAAGUGGCACCGAGUCGGUGCUUUU - 3'*The resulting sgRNA contains the desired 20 nt sequence (orange) complementary to the target DNA sequence.*

IN VITRO CRISPR/CAS9 (sgHL4)**Target DNA sequence (698 bp):**

5' -

CTTCTCCTCGACGTGCACGGACTGGACCTCGCCGAGCAGGATGATGTGGTCGCCCCGCCGGGTAGCGCTCG
 TGGACCGTGCACCTCGACGACCGCGACCGCTCCGTTCGAGCACGGTCGCTCCCCGCGCGGTACGGACGAAC
 TCCCCGCCCGCGAACTTGTCCGCGGACTTGC GCGCGAAGCGCATGGCCAGGTCCGTGTGGTCCTCGCGCA
 GCACGCTCACCGCGAACTCGCCGCAACTGTGCAACACCGGGAAGGAGTTGGCCGTACGAG**GCCAGGCAGA**
CCAGTGC|CAGTGGCGGCTCCATCGAGACGGACACGAACGAACACTGGCGGTGAAACCGTGCGGGACTCCCC
 CGCGGTTCATGGGCGGTGACGAGCGCCACCCCGGCCCGCACCCGGGCCATGGCGTCCCGGAGCATTCCCT
 GGTTCGGCTGCCATCTTCGAACTCCCTAGGCGAGGCAGGTGGGCAGCTGCCACCGTTGTAGGCGACCAGC
 TCCGCGAAGATCGCCATGGCGTCACGCACGTGCCCGCCAGCAGCTCCTCGTGGGCCCGGTGCAGGTTCCG
 CGACGAGGCGCTCCGCGTCCAGCAGCCCGGCGAAGGCGCCGAGGTCCGCCCGGCGCGCGGCCTGCAAC
 GGTGTCAGCCGCGGTTCGACCGCATCGGCGGCGAGGCGCTGGACCCAGCGCAGATAGTCCCGGTTGGCG
 TCGAT

- 3'

Cas9 protein complexed with sgRNA binds to complementary DNA sequence on the target strand after recognition of the PAM site (5'-NGG-3', in red). The Cas9 nuclease cuts 3 nt upstream of the PAM site (red line).

Cleaved sequence (285 + 413 bp):

5' -

CTTCTCCTCGACGTGCACGGACTGGACCTCGCCGAGCAGGATGATGTGGTCGCCCCGCCGGGTAGCGCTCG
 TGGACCGTGCACCTCGACGACCGCGACCGCTCCGTTCGAGCACGGTCGCTCCCCGCGCGGTACGGACGAAC
 TCCCCGCCCGCGAACTTGTCCGCGGACTTGC GCGCGAAGCGCATGGCCAGGTCCGTGTGGTCCTCGCGCA
 GCACGCTCACCGCGAACTCGCCGCAACTGTGCAACACCGGGAAGGAGTTGGCCGTACGAG**GCCAGGCAGA**
CCAGTGC - 3'

+

5' -

CAGTGGGCGGCTCCATCGAGACGGACACGAACGAACACTGGCGGTGAAACCGTGCGGGACTCCCCGCGGTC
 ATGGGCGGTGACGAGCGCCACCCCGGCCCGCACCCGGGCCATGGCGTCCCGGAGCATTCCCTGGTCGGC
 TGCCATCTTCGAACTCCCTAGGCGAGGCAGGTGGGCAGCTGCCACCGTTGTAGGCGACCAGCTCCGCGA
 AGATCGCCATGGCGTCACGCACGTGCCCGCCAGCAGCTCCTCGTGGGCCCGGTGCAGGTTCCGCGACGA
 GCGGCTCCGCGTCCAGCAGCCCGGCGAAGGCGCCGAGGTCCGCCCGGCGCGCGGCCTGCAACGGTGT
 AGCCGCGGTCGACCGCATCGGCGGCGAGGCGCTGGACCCAGCGCAGATAGTCCCGGTTGGCGTCGAT -

3'



